

Scheme: Sperm-hyaluronan assay

Distribution: Q3/4 - 2021/2022

Lab name: LabMagister

Lab number: 971

Reviewed by:


Date:

Action taken:


I Dr Mikus Endre on behalf of LabMagister have read and agreed to the terms and conditions outlined by EQASRM.

Sample code	HB393	HB394	HB395	HB396
Description	Semen sample I	Semen sample II	Semen sample III	Semen sample IV
Number of returns	4	4	4	4
Responses for normal maturity/physiological function	Yes 3/4 No 1/4	Yes 3/4 No 1/4	Yes 1/4 No 3/4	Yes 3/4 No 1/4
Your result (% sperm bound to slide)	71	76	67	74
Does this indicate normal maturity/physiological function?	No	No	No	No
Comment	Acceptable	Acceptable	Acceptable	Acceptable
<b>Overall statistic</b>				
Mean	80.0	82.8	75.0	82.8
SD	6.5	6.1	5.6	6.2
Range	[71-86]	[76-90]	[67-80]	[74-88]
Your bias	-11.3%	-8.2%	-10.7%	-10.6%
Points to note	HB393: Semen sample I after 5 minutes incubation. HB394: Sample I after 10 minutes incubation. HB395: Semen sample II after 5 minutes incubation. HB396: Semen sample II after 10minutes incubation.			


### HB393




### HB394




### HB395



### HB396



### Bias



You can alter the date range of results by grabbing and sliding the ends of this section.

